



J. N. CHEMICALS

STERISCOPE

CONCENTRATED DISINFECTANT FOR STERILIZING INSTRUMENTS

Emergency Disinfection

Non Rusting

Total Sterilisation

Economical




Rapid action

**No wastage
Sol. can be prepared
as per need**

Non Staining



**| VIRUCIDAL- HBV & HIV | FUNGICIDAL | TUBERCULOCIDAL |
| SPORICIDAL | BACTERICIDAL | PSEUDOMONACIDAL |**



STERISCOPE

THE TOTAL DISINFECTING AND STERILIZING CONCENTRATE

SUITABLE FOR ALL TYPES OF INSTRUMENTS

- Metal
- Glass
- Rubber
- Plastic objects
- Flexible fibre optic & Rigid Endoscopes
- Dental instruments.

BROAD SPECTRUM

- Sporicidal
- Virucidal
- Fungicidal
- Bactericidal
- Tuberculocidal
- Pseudomonacidal

CAN BE USED IN ALL SECTIONS OF HOSPITAL

- I.C.U.
- Surgery
- Diagnostics Centres
- Endoscopy
- Anaesthetics
- Dentistry
- Laboratory etc.

INSTRUCTION FOR USE :

To get 10% diluted solution, take 1 part of concentrated "STERISCOPE" and 9 parts of water. (i.e. 500 ml. "STERISCOPE" into 4.5 Liters of water.)

DILUTED STERISCOPE HAS 14 DAYS ACTIVE LIFE.

DISINFECTION:

Immerse precleaned instruments completely for a minimum of 10 minutes.

Remove the instruments and rinse thoroughly with sterile water.

COLD STERILISATION:

Immerse completely precleaned instruments for 5 hours. Remove the instruments and rinse thoroughly with sterile water.

COMPOSITION :

(Each 100 gms. contains)
 Glutaraldehyde 7.0 g.
 (Eq. to Glutaraldehyde USP 14%)
 1,6 Dihydroxy - 2,5 Dihexane 8.2 g.
 (Chemically bound formaldehyde)
 Polymethyl Urea Derivatives 17.6 g.
 Base aqueous Q.S.

PRECAUTIONS :

SKIN : Avoid direct contact.

EYES : In case of contact with eyes, flush with water and seek medical assistance.

FOOD : Avoid contamination.

INSTRUMENTS : Do not immerse carbon steel instruments for over 24 hours. Solution contains **rust inhibitor** but will not prevent or remove rust caused by poor plating or improper cleaning.



DIRECTION FOR USE:

Microorganisms	Exposure Time / Concentration (%)				
	5 min.	15 min.	30 min.	1 hr.	2 hrs.
Bacteris, Fungi incl. M. terrae	-	3.0%	2.0%	1.0%	-
MRSA / EHEC / Helicobacter pylori	-	3.0%	2.0%	1.0%	-
Virucidal against enveloped Viruses (incl. HBV, HIV, HCV)	0.75%	-	-	-	-
Papova Virus	-	-	4.0%	2.0%	1.0%
Rota Virus	0.75%	-	-	-	-

Always prepare the solution with cold water (max. room temperature)

Manufactured in India by:

J. N. CHEMICALS

C-1/8 G.I.D.C. industrial Estate, Phase-I,
 Vatva, Ahmedabad - 382 445, Gujarat, (India)

E-mail: jnchemicals@hotmail.com

Website: www.jnchemicals.co.in